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HOMOLOGE HORMONE:

hTSH	1 mIU/ml	WHO 68/38
hLH	500 mIU/ml	WHO 2. IRP HMG
hFSH	1000 mIU/ml	WHO 2. IRP HMB

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Hersteller:
Hubei MEIBAO Biotechnology Co., Ltd.
No.5 Hongfu Road, Shayang Economic Development Zone,
448200 Jingmen City, People's Republic of China
www.meibaobio.com

EU Importeur:
Pharma Peter GmbH,
Tarpenring 12, D-22419 Hamburg, Germany
Tel: +49 40 537 188 80

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HCG Pregnancy Test (Midstream)

medical device

(For Self-testing and In Vitro diagnostic use only)

easy quick precise

From the first day of the missed period.

1 midstream test with monoclonal antibodies

Sensitivity: 25 mIU HCG/ml urine

(A) INTENDED USE

Be Sure is a reliable pregnancy test that is easy to use and read. Be Sure HCG Pregnancy Test is a self-performing immunoassay designed for the qualitative determination of human chorionic gonadotropin (hCG) in urine for early detection of pregnancy. From the first day of your due period, it can be used at any time of the day or night. Please read the instructions for use carefully and completely before you do the test.

(B) INTENDED USER

The HCG Pregnancy Test is intended to be used at home for self-testing for women.

(C) SUMMARY AND EXPLANATION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. During normal pregnancy, hCG can be detected in urine as early as 7 days following conception, doubling every 1.3 to 2 days. At the time of the last missed menstrual period, urine hCG levels are about 100 mIU/ml with peak levels of 100,000 to 200,000 mIU/ml seen at the end of the first trimester. The presence of hCG soon after conception and its subsequent increase in concentration during early gestational growth make it an ideal marker for the early detection of pregnancy.

(D) PRINCIPLE

The HCG Pregnancy Test test is a rapid qualitative one step assay for the detection of hCG in urine. The method employs a unique combination of monoclonal dye conjugate and polyclonal-solid phase antibodies to selectively identify the hCG in the test samples with an extremely high degree of sensitivity. In less than 5 minutes, levels of hCG as low as 25 mIU/ml can be detected.

As the test sample flows through the absorbent device, the labeled antibody-dye conjugate binds to the hCG forming an antibody-antigen complex. This complex binds to the anti-hCG antibody in the positive reaction zone (“T” area) and produces a pink-purple colored band when the hCG concentration is greater than 25 mIU/ml. In the absence of hCG, there is no line in the positive reaction zone. Unbound conjugate binds to the reagents in the control zone (“C” area), producing a pink-purple band, demonstrating that the reagents are functioning properly.

(E) REAGENTS

The hCG pregnancy test stick per foil pouch.

Ingredients: Test device comprised colloidal gold coated with 1.2 mg/ml goat anti mouse, 4 mg/ml mouse anti α hCG antibody and 10 ug/ml mouse anti β hCG antibody.

(F) MATERIALS PROVIDED

Each pouch contains:

1. One HCG Pregnancy Test (Midstream)

2. Desiccant


Each box contains

1. One HCG Pregnancy Test foil pouch

2. Package insert

Additional materials required: timer or stopwatch, possibly a cup to collect the urine.

test stick

(G) STORAGE AND STABILITY

Store test strip at 4 °C–30 °C (room temperature). Avoid sunlight. The test is stable until the date imprinted on the pouch label. Do not freeze. Do not use the test after the date printed on the package. Do not use if foil pouch is damaged.

Shelf life: 36 months

KEEP OUT OF THE ACCESSIBLE FOR CHILDREN!

(H) STORAGE AND SHELF LIFE AFTER FIRST OPENING

The midstream should be used within 1 hour after the tear of the aluminum foil pouch; if the temperature is higher than 30 °C or in high humidity environment, it should be used immediately.

(I) LIMITATION OF THE PROCEDURE

- Do not use internally.
- Alcohol may interfere the test result. It is not recommended using the test after drinking.
- Occasionally specimens containing less than 25 mIU/ml for urine also yield positive results.
- A very early pregnancy containing an extremely low concentration of hCG can give a negative result. In this case, another specimen should be obtained at least 48 hours later and tested.
- hCG levels may remain detectable for several weeks after normal delivery, delivery by caesarean section, spontaneous abortion or therapeutic abortion.
- In cases where very high levels of hCG are present (>500,000 mIU/ml) a false negative result can occur due to a “Prozone” effect. If pregnancy is still suspected, simply dilute specimen 1:1 with deionized water and retest.
- If a urine sample is too dilute (ie: low specific gravity) it may not contain a representative level of hCG. If pregnancy is still suspected, a first morning urine sample should be obtained from the user in 24–48 hours and retested.
- As is true with any diagnostic procedure, the user should evaluate data obtained by the use of this kit in light of other clinical information and consult to the physicians for the final diagnosis of pregnancy before make any decision of medical relevance.

(J) ASSAY PROCEDURE

1. DETERMINATION OF TEST DATE


The test could be used from the first missing date.

2. SPECIMEN COLLECTION AND HANDLING

The Be Sure Pregnancy Test is formulated for use with fresh urine specimens. The test should be used right after the specimen collection. Urine cup should be used to collect specimens, and the urine does not require any special pretreatment. To collect the sample, you can use a urine cup or urinate directly onto the absorbent tip of the test stick. The urine does not have to be specially pre-treated. Heavy fluid intake can dilute the pregnancy hormone and falsify the test result.

3. TEST PROCEDURE

- Remove the test stick from the foil pouch
- Remove the pink cap to expose the absorbent tip
- Put the pink cap onto the end of the teststick
- Hold the stick by the thumb grip with the exposed absorbent tip pointing downward. Urinate on the suction tip for at least 10–15 seconds until it is thoroughly wet. Alternatively, you can collect the urine in a cup and hold the suction tip of the test stick in the urine for at least 10–15 seconds.
- Take the cap off the top and put it back on the suction tip. Place the stick on a clean and dry surface.
- Lay the test stick down on a flat surface with the windows on top while you wait for the test result.
- After just 1 minute, lines appear on the inspection window.
- Read the result in one to five minutes
DO NOT INTERPRET RESULT AFTER 10 MINUTES.
- Discard the test device after single use in a dustbin.

4. INTERPRETATION OF RESULTS

NEGATIVE	If only one pink-purple line appears in the control area, you can assume that you are not pregnant.
POSITIVE	If two pink-purple lines appear both in the control area and test area, you can assume that you are pregnant .
INVALID	Only line and quality control line does not appear or reaction line and quality control line didn't appear, judgement test paper is invalid. It is recommended that in this case the test be repeated.

(K) QUALITY CONTROL

Built in Quality Control Features

After addition of the sample these colored bands migrate along the membrane at the leading edge of the dye conjugate and are “removed” from the test strip completely. When the test is complete, you will see a pink-purple colored band in the “C” area of the test strip on negative samples and a pink-purple colored band in the “T” and “C” area on positive samples. The appearance of the CONTROL band indicates that the test strip is performing properly and serves as a procedural control.

(L) PERFORMANCE CHARACTERISTICS

1. ANALYTICAL SENSITIVITY & DIAGNOSTIC SENSITIVITY

HCG Pregnancy Test is an extremely rapid test used to detect the presence of hCG in urine qualitatively at levels as 25mIU/ml. The analytical sensitivity and diagnostic sensitivity are all 25mIU/ml.

2. ACCURACY

Test accuracy was assessed by comparing commercial available Syntron HCG test using 160 urine samples comprised of 90 patients with a clinical diagnosis of non-pregnant, and 70 diagnosed pregnant. HCG Pregnancy Test had 70 positive results and 90 negative results, showing 100% correlation with both Syntron HCG test and the physicians' clinical interpretation. The test accuracy is more than 99%.

3. REPRODUCIBILITY

HCG standard solution was added to normal male urine to achieve concentrations at 0 mIU/ml, 12.5 mIU/ml (50% below cutoff level), 25 mIU/ml (cutoff level), 37.5 mIU/ml (50% above cutoff level), and 50 mIU/ml (100% above cutoff level). Each specimen, at each concentration of analyte, was tested four times daily, in duplicate, for five consecutive days. A total of 40 samples at each concentration were tested. The result shows that the reproducibility of the HCG Pregnancy Test is 100%.

4. VARIABILITY

HCG standard solution was added to normal male urine to achieve concentrations at 0 mIU/ml, 12.5 mIU/ml (50% below cutoff level), 25 mIU/ml (cutoff level), 37.5 mIU/ml (50% above cutoff level), and 50 mIU/ml (100% above cutoff level). Test sample from three lots were obtained and each specimen, at each concentration of analyte, was tested four times daily, in duplicate. A total of 120 samples were tested. No variable result was observed.

5. SPECIFICITY

Menopausal Urine: A study was performed using urine specimens from 20 postmenopausal women. These specimens were chosen because urine from postmenopausal women frequently interferes with pregnancy tests due to cross reactivity with other gonadotropin hormones. All 20 urine specimens were negative when tested with the HCG Pregnancy Test.

Potentially interfering substances were added to the urine at hCG levels of 0 and 25 mIU/ml.

In any case, no interference with the HCG pregnancy test was expected:

Acetaminophen	20 mg/dl	Caffeine	20 mg/dl	Hemoglobin	20 mg/dl
Acetylsalicylic Acid	20 mg/dl	Genetic Acid	20 mg/dl	Ampicillin	20 mg/dl
Ascorbic Acid	20 mg/dl	Glucose	20 mg/dl	Tetracycline	20 mg/dl
Atropine	20 mg/dl				

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
INDEX OF SYMBOLS	PRODUCT NAME Be Sure HCG Pregnancy Test (Midstream) PZN: 00746538						
	Manufacturer		Authorized representative in the European Community/ European Union		Date of manufacture		Use-by date
	Batch code		Model number MSS-A-1		Country of manufacture „CHINA“		EU Importer
	Distributor		Do not re-use		Consult instruction for use or consult electronic instructions for use		Caution
	Temperature limitation		Keep dry		Humidity limitation		Keep away from sunlight
	Stacking layer: limit 6 layers in a carton		Medical device				

Manufacturer: Hubei MEIBAO Biotechnology Co., Ltd. No.5 Hongfu Road, Shayang Economic Development Zone, 448200 Jingmen City, People's Republic of China www.meibaobio.com	EU Importer: Pharma Peter GmbH, Tarpenning 12, D-22419 Hamburg, Germany Tel: +49 40 537 188 80
Distributor: CANEA Pharma Chemisch Pharmazeutische Vertriebsgesellschaft mbH Tarpenning 12, D-22419 Hamburg, Germany Tel: +49 40 537 188 84	Manufacturer: Hubei MEIBAO Biotechnology Co., Ltd. No.5 Hongfu Road, Shayang Economic Development Zone, 448200 Jingmen City, People's Republic of China www.meibaobio.com

Representative: Shanghai International Holding Corp. GmbH (Europe) Eiffestrasse 80, D-20537 Hamburg, Germany Tel: +49 40 251 31 75	EU Importer: Pharma Peter GmbH, Tarpenning 12, D-22419 Hamburg, Germany Tel: +49 40 537 188 80
Representative: Shanghai International Holding Corp. GmbH (Europe) Eiffestrasse 80, D-20537 Hamburg, Germany Tel: +49 40 251 31 75	Distributor: CANEA Pharma Chemisch Pharmazeutische Vertriebsgesellschaft mbH Tarpenning 12, D-22419 Hamburg, Germany Tel: +49 40 537 188 84

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

HOMOLOGOUS HORMONES:	hTSH	1 mIU/ml	WHO 68/38
	hLH	500 mIU/ml	WHO 2 nd IRP HMG
	hFSH	1000 mIU/ml	WHO 2 nd IRP HMB

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